




C-Curve PEEK



C-Curve° Ti  
Expandable

2 IMPLANTS  
1 INSTRUMENTATION



### PEEK CAGE

- LORDOSIS: 3°
- RADIOLUCENT
- ONE-STEP ANTI-BACKOUT SYSTEM



### TITANIUM CAGE

- EXPANDABLE CAGE FROM 3° TO 11°
- ONE-STEP ANTI-BACKOUT SYSTEM



### DOUBLE THREADED SCREW

- ANGULATION: 35°
- CYLINDRICAL HEAD HELPING GUIDING DURING SCREWING

#### CAGE RANGE

**Footprints** (Width X Depth)  
 15x12mm (PEEK only)  
 17x12mm (PEEK & Titanium)  
 19x14mm (PEEK & Titanium)  
 22x15mm (PEEK & Titanium)

**Heights:** 5 to 8mm  
 (by 1mm increment)

#### SCREW RANGE

Ø3.5mm X L13mm  
 Ø3.5mm X L15mm  
 Ø3.5mm X L17mm

#### STERILE PACKAGING

Cage packaged per unit  
 1 Screw box = 2 X Screws