

## PASS LP – Realignment connector


### 1 Functionality

- The realignment connector are designed in 2 sizes : 5.5 and 6.0mm
- Solid anchor point – 360° capture
- Constraint screw polyaxiality. Renders the screw monoaxial.
- Not symmetrical thus not reversible.
- 9 mm offset between the axis of the screw head and the axis of the rod.

### 2 Specific instrumentation

- None.

### 3 Illustration



### 4 References

- B02235510 : for rod Ø5,5mm
- B02236010 : for rod Ø6,0mm